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Modul BM5

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File Edit View Process Tools Help

HERMES Search: Temp-30-Days

Patient: Patient Name, Patient ID, Gender, Append Wildcard, Open Folders

Study: Description, Accession, Study ID, Append Wildcard, Open Folders

Series: Description, Body Part, Modality, Append Wildcard, Apply Detail

Date Restriction: None, Study, Series, Today, Yesterday, Last Week, Last Month, Other, From, to

Series Detail: Std Modality, All, CR, CT, DS, MR, NM, OT, PT, RT, SC, SR, US, Sub Modality, All, Print, Static, WB, Dyn, Gated, TomoAcq, Trans, Coro, Sagg, GTomoAcq, TextDml

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

Standard ROI definieren

The screenshot displays the HERMES Hybrid Viewer interface. On the left, there are three main image panels showing whole-body scans. The top panel is labeled 'Phantom_1', the middle 'Tue Apr 12 2016 15:30:40 WB1 POST1_EM2', and the right 'Wed Apr 13 2016 1 WB2 ANT2_EM2'. A red rectangle labeled 'Roi 1' is drawn on the bottom-left scan. In the center, a 'ROI/VOI' window is open, showing a list of ROIs and a context menu. The menu options include 'Rename', 'Delete', 'Create Voi from Rois', 'Insert Rois Between', 'Merge Vois', 'Time Activity Curve', 'Copy Roi to', and 'Copy Roi Mirrored to'. The 'Copy Roi to' option is highlighted, and a sub-menu is visible showing a list of studies. On the right, there is a 'HERMES' logo and a panel for 'Fri Apr 15 2016 08:51:14 WB4 POST4_EM2'. At the bottom, there are several 'Select' checkboxes and a 'Layout' section.

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ROI verschieben

The screenshot displays the HERMES software interface. At the top, a list of PET scans is shown with columns for date, time, and scan name. Below this, a grid of PET scan images is visible. A red arrow points from a text box to the 'Move ROI' option in a context menu. The text box contains the following instructions:

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The context menu includes the following options:

- Move ROI
- Move All Rois
- Increase Roi Size
- Decrease Roi Size
- Edit Roi
- Rotate Roi
- Mask Roi
- Constraint Roi
- Pixel Values Roi
- Histogram Roi
- Copy Roi
- Paste Roi
- Paste Roi Mirrored
- Paste Roi To Voi
- Delete Roi
- Move Voi
- Increase Voi Size
- Decrease Voi Size
- Mask Voi
- Pixel Values Voi
- Histogram Voi
- Copy Voi
- Paste Voi
- Delete Voi
- Export
- Tools

The bottom of the interface shows a table with columns for 'i/Voi', 'Image Number', 'Cells', 'Total', 'Mean', 'Min', 'Max', 'Median', 'Deviation', and 'Surv P'. The table lists various ROIs and VOIs for different scans.

Standard ROI maskieren

Phantom_1 - Apr 12 2016 - Hermes Hybrid Viewer PDR 2.6H

File Tools Help

HERMES

Phantom_1
kkkkkk_1
Bildgebung fuer NUK-Station
Tue Apr 12 2016 15:30:40
WB1
ANT1_EM2

Tue Apr 12 2016 15:30:40
WB1
POST1_EM2

Wed Apr 13 2016 14:53:52
WB2
ANT2_EM2

Wed Apr 13 2016 14:53:52
WB2
POST2_EM2

Thu Apr 14 2016 10:16:53
WB3
ANT3_EM2

Thu Apr 14 2016 10:16:53
WB3
POST3_EM2

Fri Apr 15 2016 08:51:14
WB4
ANT4_EM2

Fri Apr 15 2016 08:51:14
WB4
POST4_EM2

Move Roi
Move All Rois
Increase Roi Size
Decrease Roi Size
Edit Roi
Rotate Roi
Mask Roi
Constraint Roi
Pixel Values Roi
Histogram Roi
Copy Roi
Paste Roi
Paste Roi Mirrored
Paste Roi To Voi
Delete Roi
Move Voi
Increase Voi Size
Decrease Voi Size
Mask Voi
Pixel Values Voi
Histogram Voi
Copy Voi
Paste Voi
Delete Voi
Export
Tools

Apply Roi Mask
Cancel Roi Mask

Roi 1

Roi 2

Roi 3

Roi 7

Roi 8

ROI/VOI Mask

slice 1 to 1

Mask Value: 0,00

☐ Scale to Maximum Value Outside Roi/Voi

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
The screenshot displays the HERMES software interface. On the left, several patient scans are shown as silhouettes with internal organ outlines. On the right, a 'ROI/VOL' window is open, showing a table of data for various regions of interest (ROIs) across different studies. Below this, a Microsoft Excel spreadsheet titled 'Dosimetriepraktikum.xlsx' is open, showing a table with columns for 'Zeit p.i. [h]', 'ANT', 'POST', and 'geom. Mitte FIA'. Red arrows point from the 'ANT' and 'POST' columns in the Excel table to the corresponding columns in the 'ROI/VOL' window table.

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