


```
=====
Acq. Operator   : SYSTEM                      Seq. Line :    3
Acq. Instrument : G-2311                     Location  :    3
Injection Date  : 18.12.2018 15:42:36         Inj       :    1
=====
```

Different Inj Volume from Sample Entry! Actual Inj Volume : 15.000 µl

Method : C:\Chem32\1\Data\Ak_405nm\sfbnsnap 2018-12-18 15-15-19\NODAGA-NCS-Ak.M (Sequence Method)

Last changed : 17.12.2018 15:08:48 by SYSTEM

Method Info : BioSep-SEC-s3000
100% PBS mit 5 mM EDTA
1,5 ml/min

```
=====
                          Area Percent Report
=====
```

```
Sorted By      :      Signal
Multiplier     :      1.0000
Dilution       :      1.0000
Do not use Multiplier & Dilution Factor with ISTDs
```

Signal 1: DAD1 C, Sig=280,4 Ref=off

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	5.567	BV E	0.3465	170.40849	7.34864	5.3916
2	6.096	VB R	0.4910	2506.39966	71.49758	79.3002
3	8.440	BBA	0.3708	483.83795	17.16958	15.3082

Totals : 3160.64610 96.01580